



## Product Data Sheet

### **BRAIN HEART INFUSION AGAR**

**Product No.** GB-DCM-00109-1A

### **INTENDED USE**

For cultivation of fastidious pathogenic bacteria, yeasts and molds.

### **PRODUCT SUMMARY**

Brain Heart Infusion Agar is a general purpose plating medium used for the isolation, cultivation, and maintenance of a variety of fastidious and non-fastidious microorganisms. Brain Heart Infusion Agar is a modification of the medium described by ROSENOW in which the brain tissue has been replaced by brain extract and the calcium carbonate by disodium hydrogen phosphate. BHI Agar supplemented with (5 to 10%) defibrinated sheep blood is used extensively for the recovery of dimorphic fungi such as *Histoplasma capsulatum* and other pathogenic fungi such as *Coccidioides immitis*. A more selective formulation containing chloramphenicol and cycloheximide is also available that will allow the recovery of pathogenic fungi while inhibiting a wide range of bacteria and saprophytic fungi. McDonough et al. demonstrated that the temperature of incubation may alter the sensitivity of some pathogenic fungi to antibiotics; it is therefore recommended that both an antimicrobial containing medium and non-selective medium be used on primary isolates at both 25°C and 35°C.

### **Product Specifications**

<b>Ingredients</b>	<b>Gms / Ltr</b>
Calf brain infusion solids	15.000
Proteose peptone	12.500
Beef heart, infusion from	10.000
Sodium chloride	5.000
Disodium hydrogen phosphate	5.000
Dextrose	2.500
Agar	2.000

**Goslar Biotech**, 255A Barking Road East Ham, London E6 1LB, United Kingdom

Email: [info@goslarbiotech.com](mailto:info@goslarbiotech.com), Website: [www.goslarbiotech.com](http://www.goslarbiotech.com)

Microorganism	ATCC	Inoculum (CFU/ml)	Growth	Recovery	Incubation Temperature	Incubation Period
Escherichia coli	25922	50-100	Luxuriant	$\geq 70\%$	25-35°C	18-24 Hours
Streptococcus pneumoniae	49616	50-100	Luxuriant	$\geq 70\%$	25-35°C	18-24 Hours
Staphylococcus aureus	25923	50-100	Luxuriant	$\geq 70\%$	25-35°C	18-24 Hours
Shigella flexneri	12022	50-100	Luxuriant	$\geq 70\%$	25-35°C	18-24 Hours
Candida albicans	10231	10-100	Luxuriant	$\geq 70\%$	25-35°C	18-24 Hours

#### PRINCIPLE

The mixture of brain and heart infusions provides organic nitrogen, carbon, and vitamins. Dextrose is the carbohydrate source. A low concentration of dextrose is used to stimulate early growth. Sodium chloride maintains the osmotic environment. Disodium phosphate is the buffering agent in this medium. With the addition of 7% defibrinated blood the medium will support the growth of a wide range of fastidious and non-fastidious organisms, the phosphate buffer will help neutralize the acids produced from the utilization of dextrose and thus maintain viability. BHI agar can be supplemented with antibiotics, varying amounts varying amounts of Sodium chloride, Yeast extract, and serum to provide a rich medium for bacteria, yeasts and pathogenic fungi.

#### INSTRUCTION FOR USE

- Dissolve 52.0 grams in 1000 ml purified / distilled water.
- Gently heat to dissolve the medium completely.
- Sterilize by autoclaving at 15 psi (121°C) for 15 minutes.
- Cool to 45-50°C. If desired add 20 unit of Penicillin and 40µg Streptomycin/ml.
- Mix well and pour into sterile Petri plates.

#### QUALITY CONTROL SPECIFICATI

Appearance of Powder: Cream to yellow homogeneous free flowing powder.  
 Appearance of prepared medium: Basal medium: Light amber colour, clear to slightly opalescent gel. After addition of 5% v/v sterile defibrinated blood: Cherry red colour, opaque gel .  
 pH (at 25°C) : 7.4±0.2



**GOSLAR BIOTECH**  
*Your Partner in Sera*

### **STORAGE**

Dehydrated powder, hygroscopic in nature, store in a dry place, in tightly-sealed containers between 25-30°C and protect from direct sunlight. Under optimal conditions, the medium has a shelf life of 4 years. When the container is opened for the first time, note the time and date on the label space provided on the container. After the desired amount of medium has been taken out replace the cap tightly to protect from hydration.

**Product Deterioration:** Do not use if they show evidence of microbial contamination, discoloration, drying or any other signs of deterioration.

### **DISPOSAL**

After use, prepared plates, specimen/sample containers and other contaminated materials must be sterilized before discarding.

**This product is for research use only.**