



Product Data Sheet

Methanol HPLC Grade

Product No. GB-CHEM-00482
Volume. 2.5 L

TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Color Test	≤ 10 APHA
UV Absorbance 400nm	≤ 0.01
UV Absorbance 280nm	≤ 0.01
UV Absorbance 260nm	≤ 0.04
UV Absorbance 240nm	≤ 0.10
UV Absorbance 235nm	≤ 0.10
UV Absorbance 230nm	≤ 0.20
UV Absorbance 220nm	≤ 0.30
UV Absorbance 210nm	≤ 0.60
UV Absorbance 205nm	≤ 1.00
Purity (GC)	≤ 99.90 %
Residue on Evaporation	≤ 0.0005 %
Water (by Karl Fischer)	≤ 0.03 %
Fluorescence 254nm	≤ 1.0 ppb
Fluorescence 365nm	≤ 1.0 ppb

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.