

Product Data Sheet

Ethyl Acetate HPLC Grade

Product No. GB-CHEM-00484
Volume. 2.5 L

TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Titrateable Acid	Free Acid \leq 0.005 % As Acetic Acid
Purity (GC Area %)	\geq 99.70 %
Water (Coulometr.)	\leq 0.03 %
Residue (Evaporation)	\leq 0.0005 %
UV-Absorbance	\leq 0.01 AT 400NM
UV-Absorbance	\leq 0.01 AT 300NM
UV-Absorbance	\leq 0.10 AT 270NM
UV-Absorbance	\leq 0.30 AT 260NM
UV-Absorbance	\leq 0.70 AT 255 NM

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.