

## Product Data Sheet

### 2-Propanol HPLC Grade

**Product No.** GB-CHEM-00485  
**Volume.** 2.5 L

TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Color (In APHA)	≤ 10 APHA
Purity (GC Area %)	≥ 99.85 %
Water (Coulometr.)	≤ 0.05 %
Residue (Evaporation)	≤ 0.0003 %
UV-Absorbance	400nm ≤ 0.01
UV-Absorbance	260nm ≤ 0.01
UV-Absorbance	230nm ≤ 0.10
UV-Absorbance	220nm ≤ 0.20
UV-Absorbance	210nm ≤ 0.70
UV-Absorbance	205nm ≤ 0.10

**Intended Use:** For research use only.

This document has been created digitally and is valid without a signature.