

Product Data Sheet

Methyl Alcohol A.R Grade

Product No. GB-CHEM-00491
Volume. 2.5 L

TEST	Specification
Appearance (Clarity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Water (by Karl Fischer)	≤ 0.50 %
Purity (GC)	≥ 99.0 %
1H NMR Spectrum	Conforms to Structure
13C NMR Spectrum	Conforms to Structure
Deuterium NMR Spectrum	Conforms to Structure
D Enrichment	≥ 98.0 atom %
13C Enrichment	≥ 99.0 atom %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.