

Product Data Sheet

Formaldehyde Solution 37% A.R Grade

Product No. GB-CHEM-00492
Volume. 2.5 L

TEST	Specification
Assay (acidimetic)	36.5-38.0 %
Assay (iodomatic)	36.5-38.0 %
Identify	Passes test
Color	≤ 10 Hazen
Appearance of solution	Passes test
Free acide (as HCOOH)	≤ 0.025 %
Density (20°C)	1.080-1.090 g/ml
Cholaride (Cl)	≤ 0.0001 %
Sulfate (SO ₄)	≤ 0.002 %
Heavy meals (as pb)	≤ 0.0002 %
Fe (Iron)	≤ 0.0001 %
Mathanol (GC)	9.0-11.0 %
Sulfate ash	≤ 0.002%

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.