



## Product Data Sheet

### 5-Sulfosalicylic acid, sodium salt

**Product No.** GB-CHEM-00497  
**Volume.** 100G

Quality Profile	Specification
Appearance (Color)	White to Faint Yellow
Appearance (Form)	Powder or Crystals
Purity (GC Area %)	≥ 95.0 %
Carbon Content	46.75 % (Theory Monosodium Salt)
Hydrogen Content	7.19 % (Theory Monosodium Salt)
Proton NMR Spectrum	Conforms To Structure
Assay (ICP)	14.9 % (Theory Monosodium Salt)

**Intended Use:** For research use only.

This document has been created digitally and is valid without a signature.