



Product Data Sheet

Acrylamide

Product No. GB-CHEM-00499
Volume. 1KG

Quality Profile	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Crystals
Free Carboxylic Acid (as Acrylic Acid)	$\leq 0.005 \%$
Purity (HPLC)	$\geq 95.0 \%$
Solubility (Color)	Colorless
Solubility (Turbidity) 2.5g plus 10ml Water	Clear
Insoluble Matter $\leq 0.005 \%$ (Water-Insolubles)	≤ 5 ppm
Heavy Metal (as Lead)	≤ 1 ppm
Iron (Fe)	≤ 3 ppm
Magnesium (Mg)	None Detected

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.