

Product Data Sheet

AEBSF

Product No. GB-CHEM-00501
Volume. 50MG

Quality Profile	Specification
Appearance (Color)	Powder
Appearance (Form)	White
Solubility (Solvent)	Water
Solubility (Conc)	49.00 - 51.00 mg/mL
Solubility (Turbidity)	Clear
Solubility (Color)	Colorless
Initial Melting Point	≥ 183 °C
Final Melting Point	≤ 191 °C
NMR (Solvent)	D2O
Identity by NMR	Consistent
Purity (HPLC)	≥ 97.0 %
Water (by Karl Fischer)	≤ 5.00 %
Elemental Anal. (%C Anhydrous)	39.60 - 40.50 %
Elemental Anal. (%N Anhydrous)	5.40 - 6.20 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.