



Product Data Sheet

Agarose A

Product No. GB-CHEM-00502
Volume. 50G

Quality Profile	Specification
Appearance (Color)	White to off White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Light Yellow
Solubility (Turbidity) 1.5 g plus 100 mL H ₂ O	Clear to Hazy
Loss on Drying	≤ 10 %
Electroendosmosis (-MR)	0.09 - 0.13
Gel Strength	≥ 1200 g/cm ²
Gel Point	34.5 - 37.5 °C
Sulfate (SO ₄)	≤ 0.15 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.