



Product Data Sheet

Amphotericin B

Product No. GB-CHEM-00503
Volume. 1G

Quality Profile	Specification
Appearance (Turbidity)	Clear
Appearance (Form)	Solution
Solubility (Color)	Yellow
pH	9.0 - 12.0
Potency	250 - 300 µg/ml
Endotoxin Level at 10 ml/L	≤ 0.5 EU/ml
Cell Culture Test	≥ 50 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.

Goslar Biotech, 255A Barking Road East Ham, London E6 1LB, United Kingdom

Email: info@goslarbiotech.com, Website: www.goslarbiotech.com