

Product Data Sheet

Bis-acrylamide

Product No. GB-CHEM-00507
Volume. 50G

Quality Profile	Specification
Appearance (Color)	White
Appearance (Form) Water insolubles Insoluble Matter	Powder
Acrylic Acid	≤ 0.005 %
Conductivity micromho, 40% solution	≤ 0.005 %
Lead (Pb)	≤ 10
Magnesium (Mg)	≤ 5 ppm
Iron (Fe)	≤ 3 ppm
Protease Detection	≤ 2 ppm
DNase, Exonuclease Detection	None Detected
NICKase, Endonuclease Detection	None Detected
RNase Detection	None Detected
Ratio	26.8 - 31.2

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.