

Product Data Sheet

Bluo-Gal

Product No. GB-CHEM-00508
Volume. 250MG

Quality Profile	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity) 50 mg/mL, DMF	Clear
Water (by Karl Fischer)	≤ 3.0 %
Carbon	43.8 - 45.4 %
Nitrogen	3.3 - 4.0 %
Suitability Suitable for ID of Lac (+) Bacterial Colonies	Pass
Purity (TLC)	≥ 98 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.