



## Product Data Sheet

### CHAPS

**Product No.** GB-CHEM-00511  
**Volume.** 1G

Quality Profile	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity) 50 mg/ml, H <sub>2</sub> O	Clear
Water (by Karl Fischer)	≤ 8 %
Solvent Content (by GC) Dimethylformamide	≤ 0.04 %
Carbon (Dry Basis)	60.5 - 64.4 %
Nitrogen (Dry Basis)	4.0 - 5.1 %
Infrared Spectrum	Conforms to Structure
Purity (HPLC)	≥ 98 %
µmhos/cm Conductivity 5% aqueous solution	≤ 50

**Intended Use:** For research use only.

This document has been created digitally and is valid without a signature.