



## **Product Data Sheet**

### **Chitosan**

**Product No.** GB-CHEM-00512  
**Volume.** 100G

<b>Quality Profile</b>	<b>Specification</b>
Appearance (Color)	Faint Beige to Beige
Appearance (Form) Powder and/Or Chips	Powder or Flakes
Assay Deacetylation	≥ 75 %
Viscosity c = 1%, 1% Acetic Acid	20 - 300 cps

**Intended Use:** For research use only.

This document has been created digitally and is valid without a signature.