

Product Data Sheet

D-Arabitol

Product No. GB-CHEM-00514
Volume. 25G

Quality Profile	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Purity (GC)	≥ 99 %
Specific Rotation (20/D) c=0.4 in Molybdate	123 – 135°C
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity) 50mg/ml in Water	Clear

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.