



Product Data Sheet

5,5'-Dithio-bis-(2 nitrobenzoic Acid)

Product No. GB-CHEM-00518
Volume. 5G

Quality Profile	Specification
Appearance (Color)	Faint Yellow to Yellow
Appearance (Form)	Powder or Chunks
Titration with NaOH	98.5 - 101.5 %
Purity (HPLC)	> 98.5 %
Solubility (Color)	Colorless to Light Yellow
Solubility (Turbidity) 10mg/ml, EtOH	Clear to Slightly Hazy
Infrared Spectrum	Conforms to Structure

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.