



Product Data Sheet

Glycine

Product No. GB-CHEM-00519
Volume. 500G

Quality Profile	Specification
Appearance (Color)	White to off White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Very Faint Yellow
Solubility (Turbidity) 200 mg/ml, H ₂ O	Clear
Insoluble Matter	≤ 0.01 %
Chloride (Cl)	≤ 70 ppm
A ₂₈₀ UV Absorption	≤ 0.15
Heavy Metals (as Lead) by USP	≤ 20 ppm
Purity (HPLC)	≥ 99 %
Suitability	pass

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.