

## Product Data Sheet

### HEPBS

**Product No.** GB-CHEM-00520  
**Volume.** 25G

Quality Profile	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity) 25 g plus 50 mL, H <sub>2</sub> O	Clear to Slightly Hazy
Infrared Spectrum	Conforms to Structure
Water (by Karl Fischer)	≤ 0.01 %
A290 UV absorbance (33% W/W)	≤ 0.10
Heavy Metals as Pb	≤ 5 ppm
Purity (Anhydrous) Titration by NaOH	≥ 99 %

**Intended Use:** For research use only.

This document has been created digitally and is valid without a signature.