



Product Data Sheet

MOPSO, sodium salt

Product No. GB-CHEM-00527
Volume. 100G

Quality Profile	Specification
Appearance (Color)	White
Appearance (Form)	Crystalline Powder
Solubility (Color)	Colorless
Solubility (Turbidity) 20 g plus 40 ml of H ₂ O	Clear
Infrared Spectrum	Conforms to Structure
Water (by Karl Fischer)	≤ 0.5 %
A290 UV absorbance 33% W/W	≤ 0.1
Heavy Metals	≤ 5 ppm
Titration by H ₂ SO ₄ (anhydrous)	≥ 99 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.