



Product Data Sheet

MU-alpha-GAL

Product No. GB-CHEM-00528
Volume. 100MG

Quality Profile	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Very Faint Yellow
Solubility (Turbidity) 50 mg/ml, H ₂ O	Clear
Proton NMR Spectrum	Conforms to Structure
Water (by Karl Fischer)	≤ 10 %
Purity (TLC)	≥ 98 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.