

## Product Data Sheet

### N-Acetyl-L-Cysteine

**Product No.** GB-CHEM-00533  
**Volume.** 25G

Quality Profile	Specification
Appearance (Color)	White to off White
Appearance (Form)	Powder
Solubility (Colour)	Colorless
Solubility (Turbidity) 100 mg/mL, H <sub>2</sub> O	Clear
Carbon	36.0 - 38.0 %
Nitrogen	8.3 - 8.9 %
Infrared Spectrum	Conforms to Structure
% Purity (TLC)	≥ 99

**Intended Use:** For research use only.

This document has been created digitally and is valid without a signature.

