

Product Data Sheet

Netropsin

Product No. GB-CHEM-00535
Volume. 25MG

Quality Profile	Specification
Appearance (Color)	Pale Yellow to Yellow
Appearance (Form)	Powder
Solubility (Solvent)	Water
Solubility (Conc)	0.95 - 1.05 mg/ml
Solubility (Turbidity)	Clear
Solubility (Color)	Colorless to Light Yellow
Water (by Karl Fischer)	≥ 0.00 %
Elemental Anal. (%C Anhydrous)	42.45 - 43.45 %
Elemental Anal. (%N Anhydrous)	27.33 - 28.33 %
NMR (Solvent)	D2O
Identity by NMR	Consistent
NMR (Residual Solvent 1)	Acetone
NMR (% Solvent 1)	≤ 2.00 %
Purity (HPLC)	≥ 98.00 %
Purity (TLC)	≥ 98.00 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.