

Product Data Sheet

Nuclear fast red

Product No. GB-CHEM-00536
Volume. 5G

Quality Profile	Specification
Appearance (Color)	Powder
Thin Layer Chromatography	Corresponds to Standard
Solubility (Color) 0.01G IN 10ML WATER:METHANOL 1:1, HOT	Red to Very Dark Red
Wavelength (1)(UV)	550 - 560 nm
Extinction (1)(UV)	≥ 275
Wavelength (2)(UV)	516 - 526 nm
Extinction Coefficient (2)(UV) (1%/1cm)	≥ 300

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.