

Product Data Sheet

Neutral red

Product No. GB-CHEM-00537
Volume. 25G

Quality Profile	Specification
Appearance (Color)	Green to Very Dark Green and Black
Appearance (Form)	Powder or Crystals
Thin Layer Chromatography	Corresponds to Standard
Solubility (Color)	Red to Very Dark Red
Solubility (Turbidity) 0.01g in 10ml Water	Conforms to Requirements
Loss on Drying	$\leq 8.0 \%$
Wavelength (1)(UV)	535 - 545 nm
Extinction (1)(UV) (1%/1cm)	≥ 800
Solvent (UV) EtOH 50 % + CH ₃ COOH	Conforms

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.