



Product Data Sheet

Nitrilotriacetic acid

Product No. GB-CHEM-00538
Volume. 100G

Quality Profile	Specification
Appearance (Color)	White to off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Light Yellow
Solubility (Turbidity) 100 mg/mL, 1 M NaOH	Clear to Very Slightly Hazy
Infrared Spectrum	Conforms to Structure
Purity (Titration by NaOH)	> 99 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.