

## Product Data Sheet

### Phytigel

**Product No.** GB-CHEM-00541  
**Volume.** 50G

Quality Profile	Specification
Appearance (Color)	White to off White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity) 10 mg/ml, H <sub>2</sub> O	Clear to Very Hazy
Congealing Temperature 1.0% Solution	27 - 32 °C
Transmittance	≥ 80 %
Plant Cell Culture Test	Pass

**Intended Use:** For research use only.

This document has been created digitally and is valid without a signature.