



Product Data Sheet

Riboflavin

Product No. GB-CHEM-00543
Volume. 500G

Quality Profile	Specification
Appearance (Color)	Yellow to Orange
Appearance (Form)	Powder
Solubility (Color)	Yellow to Orange
Solubility (Turbidity) 10 mg/ml, 0.1M NaOH	Clear
Infrared Spectrum	Conforms to Structure
Carbon	53.2 - 55.3 %
Nitrogen	14.6 - 15.2 %
Purity (HPLC)	> 98 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.