



## Product Data Sheet

### SDS

**Product No.** GB-CHEM-00545  
**Volume.** 500G

Quality Profile	Specification
Appearance (Color)	White to off White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Very Faint Yellow
Solubility (Turbidity) 200 mg/mL, H <sub>2</sub> O	Clear
Water (by Karl Fischer)	≤ 2.0 %
Carbon	49.2 - 50.7 %
Sodium (Na)	7.0 - 9.0 %
Purity (GC)	≥ 98.5 %
Purity by Titration	≥ 98.5 %

**Intended Use:** For research use only.

This document has been created digitally and is valid without a signature.