



Product Data Sheet

**Sodium phosphate,
dibasic, anhydrous**

Product No. GB-CHEM-00552
Volume. 500G

Quality Profile	Specification
Assay (acidimetric)	100.5 %
Identity	Pass
Ca (Calcium)	0.08 ppm
Cu (Copper)	< 0.005 ppm
Fe (Iron)	0.01 ppm
K (Potassium)	1.9 ppm
Mg (Magnesium)	0.01 ppm
Chromophore Purity	99.999 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.