



Product Data Sheet

**Sodium phosphate,
monobasic, anhydrous**

Product No. GB-CHEM-00553
Volume. 500G

Quality Profile	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder or Crystals
Solubility (Color)	Colorless
Solubility (Turbidity) 100 mg/ml, H ₂ O	Clear
Loss on Drying	≤ 1.0 %
Titration with NaOH	≥ 99.0 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.