



## Product Data Sheet

**Sodium sulfate,  
anhydrate**

**Product No.** GB-CHEM-00555  
**Volume.** 500G

<b>Quality Profile</b>	<b>Specification</b>
Appearance (Color)	White to off White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity) 10 g plus 100 ml H <sub>2</sub> O	Clear
Loss on Drying	≤ 0.5 %
Purity (Titration)	≥ 99.0

**Intended Use:** For research use only.

This document has been created digitally and is valid without a signature.