



Product Data Sheet

Sodium thioglycolate

Product No. GB-CHEM-00556
Volume. 100G

Quality Profile	Specification
Appearance (Color)	White to off White
Appearance (Form)	Powder
Redox Titration Titration by Iodine	≥ 96.5 %
Solubility (Color)	Colorless to Faint Purple
Solubility (Turbidity) 200 mg/ml in Water	Clear to Slightly Hazy
Infrared Spectrum	Conforms to Structure

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.