

## Product Data Sheet

### Sorbic acid

**Product No.** GB-CHEM-00557  
**Volume.** 500G

Quality Profile	Specification
Pharmacopoea Tests	Corresponds To Requirements
Dentification A	Melting Point: 132 - 136 °C
Identification C	IR
Appearance of Solution	Clear, Colorless
Aldehydes	Max. 0.15 %, Calc. As C <sub>2</sub> H <sub>4</sub> O
Residual Solvents	Corresponds
Heavy Metals	Correspond To Requirements
Sulphated Ash	Max. 0.2 %
Water	Max. 1.0 %
Assay	99.0 - 101.0 % (Anhydrous Substance)

**Intended Use:** For research use only.

This document has been created digitally and is valid without a signature.