

Product Data Sheet

TES, free acid

Product No. GB-CHEM-00564
Volume. 50G

Quality Profile	Specification
Appearance (Form)	White
Appearance (Color)	Crystalline Powder
Solubility (Color)	colourless
Solubility (Turbidity) 25 g plus 50 ml of H ₂ O	Clear
Infrared Spectrum	Conforms to Structure
Water (by Karl Fischer)	≤ 0.2 %
A290 UV absorbance (c = 33% w/w)	≤ 0.1
Heavy Metals (as Lead)	≤ 5 ppm
Titration with NaOH	≥ 99 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.