

Product Data Sheet

Tetramethylammonium chloride (TMAC)

Product No. GB-CHEM-00565
Volume. 500G

Quality Profile	Specification
Appearance (Color)	Colorless or White
Appearance (Form)	Powder or Crystals
Titration with AgNO ₃	≥ 99.0 %
Suitability Test (HPLC)	Corresponds to Standard
¹ H NMR Spectrum	Conforms to Structure
Absorbance UV-Absorbance At 220 Nm: ≤ 0.05	≤ 0.050
Absorbance UV-Absorbance At 220 Nm: ≤ 0.04	≤ 0.040
Absorbance UV-Absorbance At 220 Nm: ≤ 0.03	≤ 0.030
Absorbance UV-Absorbance At 220 Nm: ≤ 0.02	≤ 0.020
Cyclic Voltammogram	Corresponds to Standard

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.