



Product Data Sheet

TG Buffer (TRIS-GLYCINE) 10X Solution

Product No. GB-CHEM-00566
Volume. 4L

Quality Profile	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Titration Non Aqueous Titration	1.9 - 2.0 Mol/l
Ionic Concentration	0.2 - 0.3 Mol/l
Density	1.056 - 1.066
pH Neat, 25°C	8.4 - 8.8

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.