



Product Data Sheet

Triton X-405

Product No. GB-CHEM-00574
Volume. 500ML

Quality Profile	Specification
Appearance (Clarity)	Clear to Slightly Hazy
Appearance (Color)	colourless to Light Yellow and Pink
Appearance (Form)	Liquid
Proton NMR spectrum	Conforms to Structure
Solid Content	67.0 - 73.0 %
pH 1% in Isopropanol:Water(10:6)	7.0 - 9.0

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.

Goslar Biotech, 255A Barking Road East Ham, London E6 1LB, United Kingdom
Email: info@goslarbiotech.com, Website: www.goslarbiotech.com