

Product Data Sheet

Trypan blue

Product No. GB-CHEM-00575
Volume. 10G

Quality Profile	Specification
Appearance (Color)	Dark Green to Brown or Brown-Black
Appearance (Form)	Powder or Crystals or Granules
Infrared Spectrum	Conforms to Structure
Wavelength	603 - 607 nm
Extinction Coefficient	≥ 60000
Carbon	≥ 24.6 %
Nitrogen	≥ 5.0 %
Purity (HPLC)	≥ 80.0 %
Solubility (Turbidity) Opaque	Conforms
Solubility (Color)	Dark Blue

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.