



## Product Data Sheet

### Trypsin 1:300

**Product No.** GB-CHEM-00578  
**Volume.** 50G

Quality Profile	Specification
Appearance (Color)	White to Faint Beige
Appearance (Form)	Powder
Solubility (Color)	Colorless to very Faint Yellow
Solubility (Turbidity)	Clear To Slightly Hazy
Solubility (Method)	5 Mg/MI in 1 Mm HCL
Enzymatic Activity	1200 - 1800 U/Mg

**Intended Use:** For research use only.

This document has been created digitally and is valid without a signature.