



Product Data Sheet

Uracil

Product No. GB-CHEM-00580
Volume. 100G

Quality Profile	Specification
Appearance (Color)	White to off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity) 50 mg/ml, 1 M NaOH	Clear to Slightly Hazy
EMM	7.9 - 8.2
Wavelength Lambda Max	258 - 260 N/m
Purity (HPLC)	≥ 99.0 %
¹ H NMR Spectrum	Conforms to Structure
¹³ C NMR Spectrum	Conforms to Structure

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.