



Product Data Sheet

Zinc acetate dihydrate

Product No. GB-CHEM-00585
Volume. 250G

Quality Profile	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity) 10g + 100 mL H ₂ O + 1 mL HOAc	Clear
Zinc (Zn) ID Positive for Zinc	Conforms
Titration with EDTA	≥ 98 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.