



Product Data Sheet

Zinc nitrate, hexahydrate

Product No. GB-CHEM-00587
Volume. 500G

Quality Profile	Specification
Appearance (Color)	White
Appearance (Form) Titration with EDTA 21.4 - 22.5 % Zn	Powder or Crystals or Flakes
Solubility (Color)	21.4 - 22.5 %
Solubility (Turbidity) 10% in H ₂ O	Clear to Slightly Hazy
Identity X-Ray Diffraction Conforms to Structure	Identity Corresponds

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.