

## Certificate of Analysis

### Actinomycin D

**Product No.** GB-CHEM-00589  
**Volume.** 5MG

Quality Profile	Specification
Appearance ( color)	Orange to Red
Appearance (Form)	Powder
Solubility (Solvent)	Acetonitrile
Solubility (Conc)	0,90 - 1,10 mg/mL
Solubility (Turbidity)	Clear
Purity (HPLC)	≥ 95,00 %
NMR (Solvent)	CD3OD
Identity by NMR	Consistent
Cell Culture (Suitability)	Pass

**Intended Use:** For research use only.

This document has been created digitally and is valid without a signature.