

Certificate of Analysis

Alizarin red S

Product No. GB-CHEM-00590

Quality Profile	Specification
Appearance (color)	Powder
Appearance (Form)	Yellow-Orange to Dark Orange
Appearance (Turbidity) 1 mg/mL, H ₂ O	Clear to Very Hazy
Carbon	25.0 - 49.1 %
Infrared Spectrum	Conforms to Structure
E1%	100 - 300
Wavelength	592 - 598 nm

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.