



Certificate of Analysis

Azure A chloride salt

Product No. GB-CHEM-00593
Volume. 25G

| Quality Profile | Specification |
|--|--|
| Appearance (color) | Green to Very Dark Green |
| Appearance (Form) | Powder or Crystals |
| Infrared Spectrum | Conforms to Structure |
| Solubility (Color) at 4mg Plus 4ml of Water | Dark Blue to Very Dark Blue and Dark Purple to Very Dark Purple and Blue-Purple |
| Nitrogen | 10.0 - 14.4 % |
| Dye Content | ≥ 70 % |

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.