

Certificate of Analysis

Carbon disulfide

Product No. GB-CHEM-00601
Volume. 500ML

Quality Profile	Specification
Appearance (color)	Colorless
Appearance (Form)	Liquid
Infrared Spectrum	Conforms to Structure
GC (area %)	≥ 99.90 %
Water (by Karl Fischer)	≤ 0.05 %
Residue on Evaporation	≤ 0.002 %
Color Test	≤ 10 APHA
Sulfur Dioxide	Pass
Hydrogen Sulfide	Pass
Meets ACS Requirements	Conforms

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.