

## Certificate of Analysis

### Cholesterol

**Product No.** GB-CHEM-00602  
**Volume.** 50G

Quality Profile	Specification
Appearance ( color)	White to Off White
Appearance (Form)	Powder
Infrared Spectrum	Conforms to Structure
Solubility Solubility in Alcohol Current NF	Pass
Residue on ignition (Ash) Titration with NaOH Current NF	≤ 0.1 %
Acidity	Pass
Loss on Drying	≤ 0.3 %
Specific Rotation (c = 2, 1,4-Dioxane, 25°C) Current NF	-38 - -34 °
Purity (GC)	≥ 92.5 %

**Intended Use:** For research use only.

This document has been created digitally and is valid without a signature.