

Certificate of Analysis

Glucose oxidase

Product No. GB-CHEM-00617

Quality Profile	Specification
Appearance (Color)	Yellow
Appearance (Form)	Powder
Protein content (Biuret)	≥ 60.0 %
Enzymatic activity	≥ 100,000 units/g solid
Extraneous activities	Catalase: ≤ 10 u/mg protein Galactose oxydase ≤ 3.0 % Glycogenase ≤ 0.2 % Invertase ≤ 0.5 % Maltase ≤ 0.3 % Protease ≤ 0.02 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.