

Product Data Sheet

Iron chloride (Ferric chloride), tetrahydrate

Product No. GB-CHEM-00624
Volume. 250MG

Quality Profile	Specification
Appearance (Color)	Faint Green to Green to Yellow
Appearance (Form)	Powder or Crystal or Chunk(s)
Redox Titration with KMnO ₄	97.5 - 102.5 %
Solubility (Color)	Colorless to Green to Yellow
Solubility (Turbidity) 5% in 2% HCl	Clear
ICP Major Analysis	Confirmed

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.